

Human Interleukin-2 protein (Recombinant) {Biotin} (STJP001800)

STJP001800

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Interleukin-2-protein was developed from cho cells na na. For use in research applications.
Host/Source	CHO cells NA

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 Mu m filtered PBS solution. NA
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C C. Upon reconstitution, this cytokine can be stored in
Instruction	working aliquots at °C-°C C for one month, or at -20°C C for six months, with a carrier protein without detectable

TARGET INFORMATION

Gene ID	3558
Gene Symbol	IL2
Uniprot ID	IL2_HUMAN
Immunogen	
Sequence	