

Human 2B4/CD244/SLAMF4 Protein protein (Recombinant) (STJP001802)

STJP001802

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human 2B4/CD244/SLAMF4 Protein-protein was developed from hek293 cells na na. For use in research applications.
Host/Source	HEK293 cells NA

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C C.

TARGET INFORMATION

Gene ID	51744
Gene Symbol	CD244
Uniprot ID	CD244_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081