

SARS-COV-2 Spike S1 (Omicron B.529) Protein protein (Recombinant) {Biotin} (STJP001803)

STJP001803

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-SARS-COV-2 Spike S1 (Omicron B.529) Protein-protein was developed from hek293 cells na na. For use in research applications.
Host/Source	HEK293 cells NA

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C C.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence