

HLA-A*02:01 P53 peptide complex peptide (Recombinant) (STJP001818)

STJP001818

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Recombinant-HLA-A*02:01 P53 peptide complex is synthetically produced from the HLA-A*02:01 AA25-305/B2M protein AA2-119 sequence and is suitable for use in western blot applications.
Host/Source	HEK293 cells NA

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4 NA
Purification	
Dilution Range	purity was determined to be >95% by nonreducing PAGE and HPLC (see below methods) NA
Storage	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C C. NA
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	