

Mouse CD10 protein (Recombinant) (STJP001882)

STJP001882

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse CD10-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	0.2 Mu m filtered solution in 50mM Tris, 100mM NaCL, pH8.0.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C C for six months, with a carrier protein without detectable
Instruction	loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	4311
Gene Symbol	MME
Uniprot ID	NEP_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081