

human TNFSF13B/BAFF/CD257 protein (Recombinant) (N-6His) (STJP001909)

STJP001909

GENERAL INFORMATION



Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification
Dilution Range > 92% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID [10673](#)
Gene Symbol [TNFSF13B](#)
Uniprot ID [TN13B_HUMAN](#)
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081