

Human PDGF-AA protein (Recombinant) (No-Tag) (STJP001943)

STJP001943

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	Pichia

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20 mM NaAc, pH 4.5
Purification	
Dilution Range	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5154
Gene Symbol	PDGFA
Uniprot ID	PDGFA_HUMAN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081