

Mouse IL19 protein (Recombinant) (STJP002009)
STJP002009

GENERAL INFORMATION


Product Type Proteins
Short Description
Host/Source

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification
Dilution Range > 95% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 329244
Gene Symbol IL19
Uniprot ID IL19_MOUSE
Immunogen Sequence



Recombinant Mouse IL-19 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 50-55 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081