

## Human Angiopoietin-like 4/ANGPTL4 protein (Recombinant) (N-His) (STJP002072)

STJP002072

### GENERAL INFORMATION


|                          |              |
|--------------------------|--------------|
| <b>Product Type</b>      | Proteins     |
| <b>Short Description</b> |              |
| <b>Host/Source</b>       | HEK293 cells |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Concentration</b>       |  |
| <b>Formulation</b>         | Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4.                                  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      | > 95% by SDS-PAGE.   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                           |             |
|---------------------------|-------------|
| <b>Gene ID</b>            | 51129       |
| <b>Gene Symbol</b>        | ANGPTL4     |
| <b>Uniprot ID</b>         | ANGL4_HUMAN |
| <b>Immunogen Sequence</b> |             |



---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081