

Mouse Fas ligand/APTL/CD95 ligand/CD178 protein (Recombinant) (N-His) (STJP002447)

STJP002447

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	Yeast

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	Please contact us for more information.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14103
Gene Symbol	Faslg
Uniprot ID	TNFL6_MOUSE
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081