

Human Angiopoietin-like 3/ANGPTL3 protein (Recombinant) (C-His) (STJP002462)

STJP002462

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, 0.05% CHAPS, pH 7.4.
Purification	
Dilution Range	> 90% as determined by HPLC
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27329
Gene Symbol	ANGPTL3
Uniprot ID	ANGL3_HUMAN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081