

Human ALK-4/ACVR1B protein (Recombinant) (C-hFc & his) (STJP002553)

STJP002553

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES


Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 97% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	91
Gene Symbol	ACVR1B
Uniprot ID	ACV1B_HUMAN
Immunogen Sequence	



Recombinant Human ALK-4/ACVR1B Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 46 kDa.



Immobilized Human ACVR1B at 0.5 Mu g/mL (100 Mu L/well) can bind Human ACVR2B with a linear range of 2.0-266.1 ng/mL.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081