

Human FGFR-4/CD334 protein (Recombinant) (C-hFc & His) (STJP002578)

STJP002578

GENERAL INFORMATION





Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification
Dilution Range > 95% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 2264
Gene Symbol FGFR4
Uniprot ID FGFR4_HUMAN
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081