

Human SLAMF1/CD15 protein (Recombinant) (C-hFC) (STJP002630)

GENERAL INFORMATION

| | |
|--------------------------|--------------|
| Product Type | Proteins |
| Short Description | |
| Host/Source | HEK293 cells |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | |
| Formulation | Lyophilised from a 0.2 Mu m filtered solution of PBS, pH 7.4. |
| Purification | |
| Dilution Range | > 97% by SDS-PAGE. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|-------------|
| Gene ID | 6504 |
| Gene Symbol | SLAMF1 |
| Uniprot ID | SLAF1_HUMAN |
| Immunogen Sequence | |

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081