

Human IFNA6/IFN-alpha 6/LeIF K protein (Recombinant) (C-6xHis) (STJP002634)

STJP002634

GENERAL INFORMATION

Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.2 M µm filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3443
Gene Symbol	IFNA6
Uniprot ID	IFNA6_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081