

## Human Fc-gamma RIII alpha/CD16a protein (Recombinant) (C-hFc & His) (STJP002689)

STJP002689

### GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	
<b>Host/Source</b>	HEK293 cells

### PRODUCT PROPERTIES


<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
<b>Purification</b>	
<b>Dilution Range</b>	> 95% by SDS-PAGE.
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	2214
<b>Gene Symbol</b>	FCGR3A
<b>Uniprot ID</b>	FCG3A_HUMAN
<b>Immunogen Sequence</b>	



Recombinant Human Fc-gamma RIII alpha/CD16a Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 71 kDa.



Immobilized recombinant Human FCGR3A at 1 Mu g/mL (100 Mu L/well) can bind FCGR3A Mouse mAb with a linear range of 0.13-7.8 ng/mL.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081