

Human FGL1/HP-041 protein (Recombinant) (C-His) (STJP002744)

STJP002744

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.22 Mu m filtered solution in 20mM Tris, 300mM NaCl, 10% Glycerol, pH 8.0
Purification	
Dilution Range	> 90% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2267
Gene Symbol	FGL1
Uniprot ID	FGL1_HUMAN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081