

Human HGK2/CD167a protein (Recombinant) (C-His) (STJP002764)

STJP002764

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES


Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	780
Gene Symbol	DDR1
Uniprot ID	DDR1_HUMAN
Immunogen Sequence	



Recombinant Human HGK2/CD167a Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 55-65kDa.



Immobilized recombinant Human DDR1 at 1 Mu g/mL (100 Mu L/well) can bind DDR1 Rabbit pAb with a linear range of 1-200 ng/mL.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081