

## Mouse Sonic hedgehog N-product/SHH protein (Recombinant) (C-His) (STJP002787)

STJP002787

### GENERAL INFORMATION

|                          |              |
|--------------------------|--------------|
| <b>Product Type</b>      | Proteins     |
| <b>Short Description</b> |              |
| <b>Host/Source</b>       | HEK293 cells |

### PRODUCT PROPERTIES

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| <b>Concentration</b>       |                                                                                                |
| <b>Formulation</b>         | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.                                 |
| <b>Purification</b>        |                                                                                                |
| <b>Dilution Range</b>      | > 97% by SDS-PAGE.                                                                             |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                           |           |
|---------------------------|-----------|
| <b>Gene ID</b>            | 20423     |
| <b>Gene Symbol</b>        | Shh       |
| <b>Uniprot ID</b>         | SHH_MOUSE |
| <b>Immunogen Sequence</b> |           |

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081