

Mouse IL-17RA/CD217 protein (Recombinant) (C-hFc) (STJP002810)

STJP002810

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES


Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16172
Gene Symbol	Il17ra
Uniprot ID	I17RA_MOUSE
Immunogen Sequence	



Recombinant Mouse IL-17RA/CD217 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 67-77kDa.



Immobilized Mouse IL-17A/CTLA-8 at 2 Mu g/mL (100 Mu L/well) can bind Mouse IL-17RA/CD217 with a linear range of 0.1-26 ng/mL.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081