

Human CD9 protein (Recombinant) (C-hFC) (STJP002856)

STJP002856

GENERAL INFORMATION

Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 µm filtered solution of PBS, pH 7.4.
Purification
Dilution Range > 97% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID [928](#)
Gene Symbol [CD9](#)
Uniprot ID [CD9_HUMAN](#)
Immunogen Sequence