

Mouse IL-12/IL-12A & IL-12B protein (Recombinant) (C-His (IL-12a) & No tag (IL-12b))
(STJP002890)
STJP002890

GENERAL INFORMATION

Product Type	Proteins
Short Description	
Host/Source	HEK293 cells


PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16159 16160
Gene Symbol	Il12a Il12b
Uniprot ID	IL12A_MOUSE IL12B_MOUSE

Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081