

Human Cystatin-A/CSTA protein (Recombinant) (No tag) (STJP002900)

STJP002900

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E.coli

PRODUCT PROPERTIES


Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1475
Gene Symbol	CSTA
Uniprot ID	CYTA_HUMAN
Immunogen Sequence	



Recombinant Human Cystatin-A/CSTA Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 14 kDa.



Recombinant Human OSM inhibits papain cleavage of a fluorogenic peptide substrate Z-FR-AMC. The IC50 value is <0.7 nM.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081