

Mouse IL-1Ra/IL-1F3/IL-1RN protein (Recombinant) (No-Tag) (STJP002910)

GENERAL INFORMATION


Product Type Proteins  
Short Description  
Host/Source E.coli

PRODUCT PROPERTIES

Concentration  
Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.  
Purification  
Dilution Range > 95% by SDS-PAGE.  
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 16181  
Gene Symbol Il1rn  
Uniprot ID IL1RA\_MOUSE  
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081