

Human FABP3/H-FABP protein (Recombinant) (N-His & No-Tag) (STJP002911)

STJP002911

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Purification	
Dilution Range	> 98% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2170
Gene Symbol	FABP3
Uniprot ID	FABPH_HUMAN
Immunogen Sequence	



Recombinant Human FABP3/H-FABP Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-20 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081