

Human S100-B protein (Recombinant) (N-hFc) (STJP002915)

STJP002915

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.2 Mu m filtered solution of 20mM PBS, 150mM NaCl, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6285
Gene Symbol	S100B
Uniprot ID	S100B_HUMAN
Immunogen Sequence	



Recombinant Rat IL-4 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 40-45kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081