

Human EGLN1/PHD2 protein (Recombinant) (C-Flag) (STJP002923)

STJP002923

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Hepes, 250mM NaCl, pH 7.5.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54583
Gene Symbol	EGLN1
Uniprot ID	EGLN1_HUMAN
Immunogen Sequence	



Recombinant Human EGLN1/PHD2 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 30-35kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081