

# Rat CSF-1/M-CSF protein (Recombinant) (C-His) (STJP002933)

STJP002933

## GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	
<b>Host/Source</b>	E.coli

## PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

## TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	Csf1
<b>Uniprot ID</b>	CSF1_RAT
<b>Immunogen Sequence</b>	



Recombinant Rat CSF-1/M-CSF Protein was resolved with SDS PAGE under reducing (R) and non-reducing (NR) conditions, showing single bands at 40-55 kDa and 75-110 kDa, respectively.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081