

## Mouse CCL11/Eotaxin protein (Recombinant) (No-Tag) (STJP002939)

STJP002939

### GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	
<b>Host/Source</b>	E.coli

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from 0.22 µg/ml filtered solution in PBS, 250mM NaCl, pH 7.4. Normally 8% Trehalose is added as protectant before lyophilisation.
<b>Purification</b>	
<b>Dilution Range</b>	PBS, 250mM NaCl, pH7.4
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	20292
<b>Gene Symbol</b>	Ccl11
<b>Uniprot ID</b>	CCL11_MOUSE
<b>Immunogen</b>	
<b>Sequence</b>	



Recombinant Mouse CCL11/Eotaxin Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081