

Mouse CXCL4/PF-4 protein (Recombinant) (C-His) (STJP002941)

STJP002941

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.22 µg/ml filtered solution in 20mM pb, 150mM NaCl, 1mM EDTA (pH 6.0). Normally 8% Trehalose is added as protectant before lyophilisation.
Purification	
Dilution Range	> 90% by SDS-PAGE.
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56744
Gene Symbol	PF4
Uniprot ID	PLF4_MOUSE
Immunogen	
Sequence	



Recombinant Mouse CXCL4/PF-4 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-20 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081