

Mouse PPIase A/PPIA protein (Recombinant) (No-Tag) (STJP002942)

STJP002942

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.22 μ g/ml filtered solution in PBS, 10% Glycerol, pH 7.5.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	268373
Gene Symbol	Ppia
Uniprot ID	PPIA_MOUSE
Immunogen Sequence	



Recombinant Mouse PPIase A/PPIA Protein was determined by SDS-PAGE with Coomassie Blue, showing bands at 15-25 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081