

## Human S100-A10 protein (Recombinant) (C-His) (STJP002958)

STJP002958

### GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	
<b>Host/Source</b>	HEK293 cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris 10% Glycerol, pH8.5
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	6281
<b>Gene Symbol</b>	S100A10
<b>Uniprot ID</b>	S10AA_HUMAN
<b>Immunogen Sequence</b>	



Recombinant Human S100A10 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 12.6 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081