

Human TGF- α protein (Recombinant) (No-Tag) (STJP002972)

STJP002972

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of 50mM Tris, 100mM NaCl, pH 8.0
Purification	
Dilution Range	> 92% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7039
Gene Symbol	TGFA
Uniprot ID	TGFA_HUMAN
Immunogen Sequence	



Recombinant Human TGF- α /TGFA Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081