

Mouse CCL4/MIP-1 Beta protein (Recombinant) (C-hFc) (STJP002985)

STJP002985

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Purification	
Dilution Range	> 90% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20303
Gene Symbol	Ccl4
Uniprot ID	CCL4_MOUSE
Immunogen	
Sequence	



Recombinant Mouse CCL4/MIP-1AA/ Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081