

Rat TNFRSF17/BCMA/CD269 protein (Recombinant) (N-hFc) (STJP002993)

STJP002993

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Purification	
Dilution Range	> 98% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence



Recombinant Rat TNFRSF17/BCMA/CD269 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 35-45 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081