

Human IFN- α WA/IFNA16 protein (Recombinant) (C-6His) (STJP003004)

STJP003004

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3449
Gene Symbol	IFNA16
Uniprot ID	IFN16_HUMAN
Immunogen Sequence	



Recombinant Human IFN- α WA/IFNA16 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081