

Human CCL3L1 protein (Recombinant) (N-His) (STJP003027)

STJP003027

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	414062/6349
Gene Symbol	CCL3L1.CCL3L3
Uniprot ID	CL3L1_HUMAN
Immunogen Sequence	



Recombinant Human CCL3L1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081