

Mouse IL-20 protein (Recombinant) (C-His) (STJP003035)

STJP003035

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM PB, 500mM NaCl, pH 8.0.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58181
Gene Symbol	IL20
Uniprot ID	IL20_MOUSE
Immunogen Sequence	



Recombinant Mouse IL-20 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-25 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081