

Mouse Galectin-1/LGALS1 protein (Recombinant) (C-hFC) (STJP003047)

STJP003047

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, 1mM DTT, 1mM EDTA, pH 7.4.
Purification	
Dilution Range	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16852
Gene Symbol	Lgals1
Uniprot ID	LEG1_MOUSE
Immunogen Sequence	



Recombinant Mouse Galectin-1/LGALS1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-20 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081