

Human SLAMF6/NTB-A/CD352 protein (Recombinant) (No tag) (STJP003056)

STJP003056

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E.coli

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114836
Gene Symbol	SLAMF6
Uniprot ID	SLAF6_HUMAN
Immunogen Sequence	



Recombinant Human SLAMF6/NTB-A/CD352 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081