

Human Granzyme A/CTLA-3/GZMA protein (Recombinant) (C-His) (STJP003063)

STJP003063

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	>95% as determined by SDS-PAGE
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3001
Gene Symbol	GZMA
Uniprot ID	GRAA_HUMAN
Immunogen Sequence	



Recombinant Human Granzyme A Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 30-35kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081