

Human CXCL14/BRAK protein (Recombinant) (C-His) (STJP003065)

STJP003065

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 M filtered solution of 0.5 M NaCl, 0.5 M Arginine, pH 7.4.
Purification	
Dilution Range	> 90% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9547
Gene Symbol	CXCL14
Uniprot ID	CXL14_HUMAN
Immunogen Sequence	



Recombinant Human CXCL14/BRAK Protein was determined by SDS-PAGE with Coomassie Blue, showing bands at 12 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081