

Mouse FGF-1/aFGF protein (Recombinant) (N-His (IL-23a) &No tag (IL-12b)) (STJP003074)

STJP003074

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.2 µg/ml filtered solution of 20 mM Tris-HCl, pH 7.5, 50 mM Na ₂ SO ₄ , 0.5 mM DTT, 1 mM EDTA.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14164
Gene Symbol	Fgf1
Uniprot ID	FGF1_MOUSE
Immunogen Sequence	



Recombinant Mouse FGF-1/aFGF Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-25 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081