

Streptococcus pyogenes M1 Cas9 protein (Recombinant) (N-his & hFC) (STJP003091)

STJP003091

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.22 μ g/ml filtered solution in 10mM Tris, 300mM NaCl, 0.1mM EDTA, 1mM DTT, 50% Glycerol, pH 7.4.
Purification	
Dilution Range	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence



Recombinant Streptococcus pyogenes M1 Cas9 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 110-160 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081