

Rat Leptin/LEP protein (Recombinant) (C-His) (STJP003094)

STJP003094

GENERAL INFORMATION


Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES


Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification
Dilution Range > 97% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 25608
Gene Symbol Lep
Uniprot ID LEP_RAT
Immunogen Sequence



Recombinant Rat LEP Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 16.1 kDa.



Immobilized Human LEPR (Catalog: STJP002741) at 5 Mu g/mL (100 Mu L/well) can bind Rat LEP (Catalog: STJP003094) with a linear range of 1-500 ng/mL.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081