

## Rat IL-18 protein (Recombinant) (C-His) (STJP003102)

STJP003102

### GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	
<b>Host/Source</b>	HEK293 cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
<b>Purification</b>	
<b>Dilution Range</b>	> 90% by SDS-PAGE.
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	29197
<b>Gene Symbol</b>	IL18
<b>Uniprot ID</b>	IL18_RAT
<b>Immunogen</b>	
<b>Sequence</b>	



Recombinant Rat IL-18 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 20-25 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081