

Human Collagen II alpha 1/COL2A1 protein (Recombinant) (STJP003110)

STJP003110

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	Human

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification	
Dilution Range	> 95% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1280
Gene Symbol	COL2A1
Uniprot ID	CO2A1_HUMAN
Immunogen Sequence	



Recombinant Human Collagen II alpha 1/COL2A1 Protein was resolved with SDS PAGE under reducing (R) and non-reducing (NR) conditions, showing single bands at 35-40 kDa and 105-120 kDa, respectively.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081