

Mouse Anti-HSA HSA antibody {Unconjugated} (STJS000003)

STJS000003

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Secondary antibodies |
| Short Description | Unconjugated mouse monoclonal anti-HSA HSA secondary antibody. For use in most research applications. |
| Applications | WB/IHC/IF/ELISA |
| Host/Source | Mouse |
| Reactivity | HSA |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | |
| Concentration | 2 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range | WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000 |
| Formulation | Liquid in PBS containing 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| |
|---------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen Sequence |