

Rabbit Anti-HSA HSA antibody {FITC} (STJS000065)

STJS000065

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Secondary antibodies |
| Short Description | FITC-conjugated rabbit polyclonal anti-HSA HSA secondary antibody. For use in most research applications. |
| Applications | WB/IHC/IF/ELISA |
| Host/Source | Rabbit |
| Reactivity | HSA |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 2 mg/mL |
| Conjugation | FITC |
| Purification | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range | IF 1:250-1:2000 FCM 1:250-1:2000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| |
|---------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen Sequence |