

Rabbit Anti-FITC FITC antibody {Biotin} (STJS000082)

STJS000082

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated rabbit polyclonal anti-FITC FITC secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	FITC

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Biotin
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	IF 1:250-1:2000 FCM 1:250-1:2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence